

Reusable bed pad

Polyester



- Partial mattress overlay
- Waterproof and Breathable
- Absorbent
- Material suitable for allergy sufferers
- Resistant to dust mites
- Virusproof
- Blood and urine resistant
- Can be disinfected with spray
- Hardly any slippage, thanks to brushed reverse
- Hardly any wrinkling
- Absorption volume: 2,674 ml/m²

size (cm)	80 x 90
item no.	3111 004 087



- Partial mattress overlay
- Waterproof and Breathable
- Absorbent
- Material suitable for allergy sufferers
- Resistant to dust mites
- Virusproof
- Blood and urine resistant
- Can be disinfected with spray
- Hardly any slippage, thanks to brushed reverse
- Hardly any wrinkling
- Absorption volume: 2,674 ml/m²

size (cm)	80 x 180 (80 x 90)
item no.	3112 004 084

Reusable bed pad

Cotton



- Partial mattress overlay
- Waterproof and Breathable
- Absorbent
- Material suitable for allergy sufferers
- Resistant to dust mites
- Anti-allergenic
- Virusproof
- Blood and urine resistant
- Can be disinfected with spray
- Can be washed up to 300 x
- Absorption volume: 3,449 ml/m²

size (cm)	85 x 90
item no.	3103 001 085

suprima GmbH
Kulmbacher Str. 31-35
95460 Bad Berneck


suprima
www.suprima-gmbh.de

Telefon 09273 9204 0
Telefax 09273 9204 55
info@suprima-gmbh.de

Bed Protection

you can rely on

suprima



Simply a better quality of life.

For a good night's sleep

Protection and comfort

Bed protection pads and fitted sheets are used to protect the mattress from soaking.

For temporary use or short-term care, simple pads made of PVC can be used.

Coated bed pads (terry or molleton) also have an absorbent function and increase lying comfort. They are placed on top of the sheets, so that in the event of micturition (wetting) only the bed pad must be washed.

In heavy, long-term and intensive care, pads made of PU (polyurethane) are preferred. They are more durable and also breathable. Pads made of PU reduce the risk of bedsores due to their elasticity. Multiple bed pads are highly absorbent and also prevent bedsores.

If the patient is motor-agitated, underlays should be fixed by side parts or bed protection with circumferential rubber should be used.

Fitted sheet

PU



- Material suitable for allergy sufferers
- Anti-allergenic
- Waterproof
- Resistant to dust mites
- Virusproof
- Blood and urine resistant
- Breathable
- Can be disinfected with spray

size (cm)	100x200 x20	140 x 200 x 20
item no.	3527 001 202	3527 001 201

Fitted sheet

PVC



- Material suitable for allergy sufferers
- Anti-allergenic
- Waterproof
- Corners heat-sealed with rubber
- Resistant to dust mites
- Virusproof
- Blood and urine resistant
- Can be disinfected with spray

size (cm)	100 x 200 x 20
item no.	3066 001 202

Draw sheet

Terry cloth



STANDARD 100

- Material suitable for allergy sufferers
- Anti-allergenic
- Waterproof
- PU coated
- Resistant to dust mites
- Germproof
- Breathable
- Elasticated
- Absorbent
- Securely fixed using side elements

size (cm)	90 x 150
item no.	3051 001 090

Terry cloth /PU coated



STANDARD 100

- Material suitable for allergy sufferers
- Anti-allergenic
- Resistant to dust mites
- Waterproof
- PU coated
- Germproof
- Absorbent
- Breathable
- Can be washed at least 100 x

sizes (cm)	100x200 x24	90 x 200 x 12
item no.	3071 001 104	3071 001 103



STANDARD 100



- Material suitable for allergy sufferers
- Anti-allergenic
- Waterproof
- PVC coated
- Absorbent
- Blood and urine resistant
- suitable for cots

sizes (cm)	100x200x20	140x200x20	70x140x15
item no.	3067 001 202	3067 001 201	3067 001 072